



# Protective suit

<b>Description</b>	Disposable protective suit with hood for use in the medical field.			
<b>Certification</b>	Meets the protection requirements of: DIN EN 14126:2003/AC:2004 DIN EN 14605:2005+A1:2009 CE-certified			
<b>Sterility</b>	Not sterilized			
<b>Storage</b>	Store in a cool, dry place and with sufficient ventilation. Protect from direct sunlight and keep away from ignition sources.			
<b>Available sizes</b>	Unit: cm			
	<b>Size</b>	<b>160-165</b>	<b>170-175</b>	<b>180-185</b>
	Total length	186	196	206
	Chest measurement	119	127	135
	Arm length	86	89	92
	Sleeve opening	18	18	18
	Foot opening	24	24	24

## Safety instructions

1. For single use only.
2. Do not use if the packaging is opened or damaged.
3. Dispose of after use in accordance with local laws and regulations.

We are looking forward to your order by e-mail: [sampling@inveox.com](mailto:sampling@inveox.com)

Due to high demand and the need to integrate and optimize various supply chains, the delivered product may vary in appearance and packaging from item in this description.

inveox GmbH • Lichtenbergstraße 8, 85748 Garching bei München, Deutschland  
www.inveox.com • Telephone: +49 (0) 89 / 57 84 76 01 • Fax: +49 (0) 89 / 41 41 60 543