


Disposable Nitrile Glove




Description	The powder-free disposable nitrile gloves are available in six different sizes and are suitable for the following uses: medical checking, mechanical maintenance, cleaning, tattooing and food preparation.
Material	Nitrile
Quantity	100 units per box, 10 boxes per carton
Carton size	34 × 25 × 25 cm (13.5 x 10 x 10 in)
Weight	3.5 g (.123 oz) per glove
Storage	Store in a cool dry place away from heat sources or sunlight

Expiration	Check packaging	
Quality standard	CE, EN455	
Available sizes	<p>6/XS 2.5–2.9 in</p> <p>7/S 2.9–3.4 in</p> <p>8/M 3.4–3.8 in</p> <p>9/L 3.8–4.2 in</p> <p>10/XL 4.2–4.6 in</p> <p>11/XXL 4.6–5 in</p>	 <p>Hand width</p>

Product characteristics

Superior Resistance

- Anti-puncture
- Hard to break



High Elasticity

- Strong adhesion
- Good flexibility



Fingertip texture

- Non-slip, increased friction
- Firm grip



To place an order, please get in touch via email: sampling@inveox.com

Due to high demand and the need to integrate and optimize various supply chains, the delivered product may vary in appearance and packaging from item in this description.

inveox GmbH • Lichtenbergstraße 8, 85748 Garching bei München, Deutschland
www.inveox.com • Telefon: +49 (0) 89 / 57 84 76 01 • Fax: +49 (0) 89 / 41 41 60 543